

Email: sales@laboratoryinstrumentindia.com

**Product Name :**Jet Attachment Investigation

Product Code:



## **Description:**

The shape of the duct may be changed by tilting and sliding moveable elements to allow students to investigate the Coanda effect of jet attachment to a wall.

A rectangular slit directs an air jet towards a Y shaped duct with two outlet passages.

By blowing air from one side or another at the jet, the airflow can be directed down either of the Y shaped passages as in a pneumatic flip-flop.

## **Technical Specification:**

Jet Attachment Investigation



**Laboratory Instrument India**