L	Email : sales@laboratoryinstrumentindia.com

Product Code :

Magnetic Damping Accessory To Add Friction On Cart	LBNY-0004-780013	
www.laboretoytintrument/hdlk.com		
Description :		
Magnetic Damping Accessory To Add Friction On Cart		
Technical Specification :		
Mechanical bench accessory 11: The Magnetic Damping Accessory can be used with either a l Accessory is attracted to the magnets built into the cart (or the bracket easily slides up and down to adjust the amount of dar The magnetic drag is caused by the induced currents in the a	e Magnetic Bumper on the Smart Cart), and the mping.	

The magnetic drag is caused by the induced currents in the aluminum track, and is proportional to the speed of the cart. This accessory will not work with plastic tracks!

Magnetic damping accessory to add friction on cart, supplied with: Magnets

Keeper

Adequate bracket.

Product Name :

İj

Laboratory Instrument India