

Email: sales@laboratoryinstrumentindia.com

Product Name:
Alcohol Burner

Product Code:



Description:

Alcohol Burner

Technical Specification:

Alcohol burner Material: Glass Capacity: 125 ml Plastic Cap

Alcohol Burner Lamps are used in scientific experiments requiring low heat sources. You can cut it into any length by the overlong wick.

Top up the burner to half to two-thirds. Use a funnel of the proper size to fill it in the fumigation room away from the open fire. Only use a match or barbecue lighter to light the burner in this position.

Do not attempt to ignite the burner from another ignited burner or take the ignited burner to another location. Never refuel a burning alcohol lamp. Use burner cover to extinguish burner flame.



Laboratory Instrument India