



Email : sales@laboratoryinstrumentindia.com

Product Name :
Voltage Regulator

Product Code :
LBNY-0002-670007



Description :

Voltage Regulator

Technical Specification :

Heavy-duty continuous use transformer.

On/Off switch with indicator lamp, integrated automatic protection circuit

Analog meter on front of unit measures exact input voltage.

Voltage: 0-400V.

Wattage: 3000W.

Frequency: 50Hz.

Voltage Regulator, Voltage Regulator Equipments, Voltage Regulator Tools, Voltage Regulator Tool Kits, Voltage Regulator Manufacturers, Voltage Regulator Suppliers from India, China, Kenya



Laboratory Instrument India