



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Voltage Regulator

**Product Code :**  
LBNY-0002-670007



**Description :**

Voltage Regulator

**Technical Specification :**

Heavy-duty continuous use transformer.  
On/Off switch with indicator lamp, integrated automatic protection circuit  
Analog meter on front of unit measures exact input voltage.  
Voltage: 0-400V.  
Wattage: 3000W.  
Frequency: 50Hz.

Voltage Regulator, Voltage Regulator Equipments, Voltage Regulator Tools, Voltage Regulator Tool Kits,  
Voltage Regulator Manufacturers, Voltage Regulator Suppliers from India, China, Kenya



**Laboratory Instrument India**