



Email : sales@laboratoryinstrumentindia.com

Product Name :
Impact Test 300Nm

Product Code :
LBNY-0005-10100013



Description :

In the field of industrial quality control, the impact test is a widely used test method with which to quickly and easily determine characteristics for a material or component analysis.

Features:

Charpy notched-bar impact test with increased work capacity up to 300Nm.
Safe experiments thanks to two-handed release of the hammer and optional protective cage.
Pendulum impact tester based on industrial standards.

Learning Objectives And Experiments:

Influence of notch shape, material and specimen temperature on the notched-bar impact work.
Determine the notched-bar impact work
Plot a notched-bar impact worktemperature diagram
Determine the notched-bar impact strength
Analyse the fracture surface characteristics.

Technical Specification :

Impact Test 300Nm

Laboratory Instrument India

