



Email : sales@laboratoryinstrumentindia.com

Product Name :
Digital Storage Oscilloscope (100MHz, 2 channel)

Product Code :
LBNY-0002-620031



Description :

Digital Storage Oscilloscope (100MHz, 2 channel)

Technical Specification :

Independent vertical scale & position control knobs for each channel.
Edge, Pulse Width, video, Slope, Alternate trigger mode.
Math functions including Add, Subtract, Multiply, Divide & 1024 point FFT
32 parameters of automatic measurements
Unique Digital Filter & waveform recorder function.
Advanced cursor modes: Manual, Auto & Track
Waveform Intensity & Grid Brightness can be adjusted
PASS/ FAIL detection, PASS/FAIL output
Built-in 100 MHz hardware frequency counter
AUTO's strategy presetable
Automatic measurement of 34 waveform
Automatic cursor tracking measurement
Low noise floor, wide vertical range:
7 inches TFT LCD, high resolution display.

Digital Storage Oscilloscope (100MHz, 2 channel), Digital Storage Oscilloscope (100MHz, 2 channel) Equipments, Digital Storage Oscilloscope (100MHz, 2 channel) Tools, Digital Storage Oscilloscope (100MHz, 2 channel) Tool Kits, Digital Storage Oscilloscope (100MHz, 2 channel) Manufacturers, Digital Storage Oscilloscope (100MHz, 2 channel) Suppliers from India, China, Kenya



Laboratory Instrument India