



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Function Generator (40MHz, 2 channel)

**Product Code :**  
LBNY-0002-620032



#### **Description :**

Function Generator (40MHz, 2 channel)

#### **Technical Specification :**

High precision, stable, pure and low distortion signals; provide high frequency square wave with quick rise and fall edges.

Users can improve their working efficiency by using the multifunctional instrument of easy operated interface, competitive specifications and humanized graphical display.

Designing, testing, troubleshooting, and repairing of electronic and electroacoustic devices, such as circuit and equipment testing in cellular and wireless systems, broadcasting, television and radio systems; semiconductor component testing etc.

Output wave: Sine, Square, Ramp, Burst, Noise, DC Arbitrary, Harmonic, Expression

Wave length: 8pts~16Mpts.

Function Generator (40MHz, 2 channel), Function Generator (40MHz, 2 channel) Equipments, Function Generator (40MHz, 2 channel) Tools, Function Generator (40MHz, 2 channel) Tool Kits, Function Generator (40MHz, 2 channel) Manufacturers, Function Generator (40MHz, 2 channel) Suppliers from India, China, Kenya

---

# Laboratory Instrument India

