



Email : sales@laboratoryinstrumentindia.com

Product Name :
Function Generator

Product Code :
LBNY-0002-630002



Description :

Function Generator

Technical Specification :

The models offer 50 built-in standard function and arbitrary waveforms with continuous, modulation, sweep and burst modes to cover almost all test requirements in basic education laboratories.

A built-in 200MHz counter with 6-digit resolution offers an easy and precise way of performing frequency, period, pulse width and duty cycle measurements.

Complete dual channels function/arbitrary waveform generator

Channel independence, coupling, track working mode

200 MSa/s sampling rate and 14-bit vertical resolution per channel

Output of 6 standard waveforms, built-in 50 kinds of arbitrary waveform

Multi modulation function: AM, DSSC - AM,FM,PM, PWM,FSK,ASK,-BPSK and logarithm/linear sweep.

Function Generator, Function Generator Equipments, Function Generator Tools, Function Generator Tool Kits, Function Generator Manufacturers, Function Generator Suppliers from India, China, Kenya



Laboratory Instrument India

