



Email : sales@laboratoryinstrumentindia.com

Product Name :
Bending Device Assembly Kit

Product Code :
LBNY-0005-10100039



Description :

They are grouped into individual assembly sequences. The base plate is covered with a transparent cover panel on which a graphical representation of the assembly structure is printed. The kit contains all the parts required to assemble a functional bending device. The parts are clearly laid out on a base plate. The surfaces of the finished parts have a gunmetal finish to prevent corrosion. Multiple trays are stackable, providing for space-saving storage. The kit is ideally deployed for student exercises during teaching, with two/three students able to work with each kit. The graphic symbolises standard and production parts, and also represents fixed and moving connections differently. The comprehensive and clearly structured instructional material is laid-out in line with modern principles and is of great help with lesson preparation. All parts are precision-manufactured on CNC machines and feature standard engineering tolerances and surface finishes.

Technical Specification :

Assembly set with components of a functional bending device (eccentric operation)
All components made from steel, precision-manufactured, with gunmetal surface finish
Instructional material incorporates action-oriented and interdisciplinary forms of teaching
Part of the technical drawing course
Assembly schematic on transparent cover panel
Multiple trays stackable
Dimensions and Weight of Bending Device Assembly Kit
l x w x h: 540 x 350 x 75 mm (tray)
Weight: 6 kg.



Laboratory Instrument India