



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Thread Types Kit

**Product Code :**  
LBNY-0005-10100042



**Description :**

This instructional kit is used for demonstration and information purposes. The thread flanks are made visible by cut-outs. A thread gauge enables the thread type and size to be determined. The kit is clearly laid out on a plastic tray. No provision is made for conducting exercises or experiments. Various bolt and nut threads are included. The well-structured instructional material enhances the informational value of the kit.

**Technical Specification :**

Thread flanks visible with cut-outs  
Thread gauge for determining thread type  
Thread type demonstration kit  
8 thread types, clearly laid out  
All parts clearly laid out on a plastic tray  
Galvanized parts  
Multiple trays stackable.

**Technical Data:**

Thread size: 24mm diameter  
Thread gauge for male and female threads  
Whitworth pipe thread  
Whitworth thread  
Dimensions and Weight of Thread Types Kit  
l x w x h: 500x350x110mm (tray)  
Weight: approx. 5kg.



**Laboratory Instrument India**