



Email : sales@laboratoryinstrumentindia.com

Product Name :
Thread Types Kit

Product Code :
LBNY-0005-10100042



Description :

This instructional kit is used for demonstration and information purposes. The thread flanks are made visible by cut-outs. A thread gauge enables the thread type and size to be determined. The kit is clearly laid out on a plastic tray. No provision is made for conducting exercises or experiments. Various bolt and nut threads are included. The well-structured instructional material enhances the informational value of the kit.

Technical Specification :

Thread flanks visible with cut-outs
Thread gauge for determining thread type
Thread type demonstration kit
8 thread types, clearly laid out
All parts clearly laid out on a plastic tray
Galvanized parts
Multiple trays stackable.

Technical Data:

Thread size: 24mm diameter
Thread gauge for male and female threads
Whitworth pipe thread
Whitworth thread
Dimensions and Weight of Thread Types Kit
l x w x h: 500x350x110mm (tray)
Weight: approx. 5kg.



Laboratory Instrument India