



Email : sales@laboratoryinstrumentindia.com

Product Name :
Oscilloscope

Product Code :
LBNY-0002-630003



Description :

Oscilloscope

Technical Specification :

2 channels, low noise floor, wide vertical range: 1mV/div~20V/div

100 MHz

AUTO's strategy presettable

Automatic measurement of 34 waveform

Automatic cursor tracking measurement

Low noise floor, wide vertical range:

Built in FFT and digital filter

System software upgrade via USB drive

7 inches TFT LCD, high resolution display.

Oscilloscope, Oscilloscope Equipments, Oscilloscope Tools, Oscilloscope Tool Kits, Oscilloscope Manufacturers, Oscilloscope Suppliers from India, China, Kenya



Laboratory Instrument India

