



Email : sales@laboratoryinstrumentindia.com

Product Name :
Oscilloscope

Product Code :
LBNY-0002-630003



Description :

Oscilloscope

Technical Specification :

2 channels, low noise floor, wide vertical range: 1mV/div~20V/div
100 MHz
AUTO's strategy presettable
Automatic measurement of 34 waveform
Automatic cursor tracking measurement
Low noise floor, wide vertical range:
Built in FFT and digital filter
System software upgrade via USB drive
7 inches TFT LCD, high resolution display.

Oscilloscope, Oscilloscope Equipments, Oscilloscope Tools, Oscilloscope Tool Kits, Oscilloscope Manufacturers, Oscilloscope Suppliers from India, China, Kenya



Laboratory Instrument India

