




Email : sales@laboratoryinstrumentindia.com

Product Name :
Cutaway Model Control Gear

Product Code :
LBNY-0005-10100048



Description :

The mechanical functionality of all components is retained in full. Each of the cutaway models is securely mounted on a base plate to which handles are attached for carrying purposes. All mechanisms are hand-driven. The range consists of a series of different cutaway models of drive system components. All components are industrially relevant and have been prepared such that individual active elements are clearly visible to the user.

Technical Specification :

Industrially-relevant part, fully functional cutaway model
Solid metal base plate, with handles
Hand-operated tabletop cutaway model for the demonstration of the functionality of a mechanical control gear
Dimensions and Weight:
l x w x h : 350 x 300 x 160 mm
Weight: approx. 10.5 kg.



Laboratory Instrument India