




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Cutaway Model Control Gear

**Product Code :**  
LBNY-0005-10100048



**Description :**

The mechanical functionality of all components is retained in full. Each of the cutaway models is securely mounted on a base plate to which handles are attached for carrying purposes. All mechanisms are hand-driven. The range consists of a series of different cutaway models of drive system components. All components are industrially relevant and have been prepared such that individual active elements are clearly visible to the user.

**Technical Specification :**

Industrially-relevant part, fully functional cutaway model  
Solid metal base plate, with handles  
Hand-operated tabletop cutaway model for the demonstration of the functionality of a mechanical control gear  
Dimensions and Weight:  
l x w x h : 350 x 300 x 160 mm  
Weight: approx. 10.5 kg.



**Laboratory Instrument India**