



Email : sales@laboratoryinstrumentindia.com

Product Name :
Nuts And Bolts Kit

Product Code :
LBNY-0005-10100049



Description :

The industrially-standard items are clearly laid out, screw-fitted into an aluminium panel. The symbols on the panel show the correct graphical representation as well as the DIN and standard designations of the items concerned. This instructional kit is used for demonstration and information purposes. No provision is made for conducting exercises or experiments. The kit is clearly laid out on a plastic tray. The well-structured instructional material enhances the informational value of the kit.

Technical Specification :

Screen-printing shows graphical representation and DIN designation
All parts clearly laid out on a plastic tray
Multiple trays stackable
Nuts and bolts demonstration kit
42 standard parts, clearly laid out on an aluminium panel
Screen-printed panel

Technical Data:

Aluminium panel, l x w: 350x255mm
Nuts: M6, M8, M10
Bolts: C4.8, K4.8, M6, M8, M10
Dimensions and Weight:
l x w x h: 500x350x110mm (tray)
Weight: approx. 5kg.



Laboratory Instrument India