



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Nuts And Bolts Kit

**Product Code :**  
LBNY-0005-10100049



**Description :**

The industrially-standard items are clearly laid out, screw-fitted into an aluminium panel. The symbols on the panel show the correct graphical representation as well as the DIN and standard designations of the items concerned. This instructional kit is used for demonstration and information purposes. No provision is made for conducting exercises or experiments. The kit is clearly laid out on a plastic tray. The well-structured instructional material enhances the informational value of the kit.

**Technical Specification :**

Screen-printing shows graphical representation and DIN designation  
All parts clearly laid out on a plastic tray  
Multiple trays stackable  
Nuts and bolts demonstration kit  
42 standard parts, clearly laid out on an aluminium panel  
Screen-printed panel

**Technical Data:**

Aluminium panel, l x w: 350x255mm  
Nuts: M6, M8, M10  
Bolts: C4.8, K4.8, M6, M8, M10  
Dimensions and Weight:  
l x w x h: 500x350x110mm (tray)  
Weight: approx. 5kg.



**Laboratory Instrument India**