




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Lever Shears Assembly Kit

**Product Code :**  
LBNY-0005-10100054



#### **Description :**

The kit contains all the parts required to assemble a functional lever shears. The graphic symbolises standard and production parts, and represents fixed and moving connections differently. All parts are precision-manufactured on CNC machines and feature standard engineering tolerances and surface finishes. The parts are clearly laid out on a base plate. They are grouped into individual assembly sequences. The base plate is covered with a transparent cover panel on which a graphical representation of the assembly structure is printed. Multiple trays are stackable, providing for space-saving storage. The surfaces of the finished parts have a gunmetal finish to prevent corrosion.

#### **Technical Specification :**

Instructional material incorporates action-oriented and interdisciplinary forms of teaching  
All lever shears components made from steel, precision-manufactured, gunmetal surface finish  
Assembly schematic on transparent cover panel  
Part of the technical drawing course  
Assembly kit with components of a functional lever shears  
Multiple trays stackable.

#### **Dimensions and Weight:**

l x w x h: 540 x 350 x 70 mm (tray)

Weight: approx. 7 kg.



**Laboratory Instrument India**