




Email : sales@laboratoryinstrumentindia.com

Product Name :
Lever Shears Assembly Kit

Product Code :
LBNY-0005-10100054



Description :

The kit contains all the parts required to assemble a functional lever shears. The graphic symbolises standard and production parts, and represents fixed and moving connections differently. All parts are precision-manufactured on CNC machines and feature standard engineering tolerances and surface finishes. The parts are clearly laid out on a base plate. They are grouped into individual assembly sequences. The base plate is covered with a transparent cover panel on which a graphical representation of the assembly structure is printed. Multiple trays are stackable, providing for space-saving storage. The surfaces of the finished parts have a gunmetal finish to prevent corrosion.

Technical Specification :

Instructional material incorporates action-oriented and interdisciplinary forms of teaching
All lever shears components made from steel, precision-manufactured, gunmetal surface finish
Assembly schematic on transparent cover panel
Part of the technical drawing course
Assembly kit with components of a functional lever shears
Multiple trays stackable.

Dimensions and Weight:

l x w x h: 540 x 350 x 70 mm (tray)

Weight: approx. 7 kg.



Laboratory Instrument India