



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Screw-Locking Devices Kit

**Product Code :**  
LBNY-0005-10100061



**Description :**

Screw-locking devices are presented in the as-fitted position on an aluminium panel. A transparent multi-compartment storage tray holds large numbers of screw-locking devices. This instructional kit is used for demonstration and information purposes. The well-structured instructional material enhances the informational value of the kit. No provision is made for conducting exercises or experiments. The kit is clearly laid out on a plastic tray.

**Technical Specification :**

Transparent storage tray with 18 compartments each holding 10 of the different screw-locking devices  
Screw-locking devices demonstration kit  
18 fitted screw-locking devices, clearly laid out on an aluminium panel  
Multiple trays stackable  
All parts clearly laid out on a plastic tray l x w x h: 500x350x110mm (tray)  
Weight: approx. 5kg

**Technical data:**

Locking devices for m6  
Aluminium panel, l x w: 238x100mm  
Dimensions And Weight Of Screw-Locking Devices Kit  
L x w x h: 500x350x110mm (Tray)  
Weight: Approx. 5kg.



**Laboratory Instrument India**