



Email : sales@laboratoryinstrumentindia.com

Product Name :
Screw-Locking Devices Kit

Product Code :
LBNY-0005-10100061



Description :

Screw-locking devices are presented in the as-fitted position on an aluminium panel. A transparent multi-compartment storage tray holds large numbers of screw-locking devices. This instructional kit is used for demonstration and information purposes. The well-structured instructional material enhances the informational value of the kit. No provision is made for conducting exercises or experiments. The kit is clearly laid out on a plastic tray.

Technical Specification :

Transparent storage tray with 18 compartments each holding 10 of the different screw-locking devices
Screw-locking devices demonstration kit
18 fitted screw-locking devices, clearly laid out on an aluminium panel
Multiple trays stackable
All parts clearly laid out on a plastic tray l x w x h: 500x350x110mm (tray)
Weight: approx. 5kg

Technical data:

Locking devices for m6
Aluminium panel, l x w: 238x100mm
Dimensions And Weight Of Screw-Locking Devices Kit
L x w x h: 500x350x110mm (Tray)
Weight: Approx. 5kg.



Laboratory Instrument India