



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Cutaway Model Variable Speed Belt Drive

**Product Code :**  
LBNY-0005-10100062



**Description :**

The mechanical functionality of all components is retained in full. Each of the cutaway models is securely mounted on a base plate to which handles are attached for carrying purposes. The range consists of a series of different cutaway models of drive system components. All components are industrially relevant and have been prepared such that individual active elements are clearly visible to the user. All mechanisms are hand-driven.

**Technical Specification :**

Industrially-relevant part, fully functional cutaway model  
Solid metal base plate, with handles  
Hand-operated benchtop cutaway model for the demonstration of the functionality of a variable speed belt drive  
Dimensions and weight:  
L x w x h: 350 x 300 x 200 mm  
Weight: approx. 7.5 kg.



**Laboratory Instrument India**