



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cutaway Model Variable Speed Belt Drive

Product Code :
LBNY-0005-10100062



Description :

The mechanical functionality of all components is retained in full. Each of the cutaway models is securely mounted on a base plate to which handles are attached for carrying purposes. The range consists of a series of different cutaway models of drive system components. All components are industrially relevant and have been prepared such that individual active elements are clearly visible to the user. All mechanisms are hand-driven.

Technical Specification :

Industrially-relevant part, fully functional cutaway model
Solid metal base plate, with handles
Hand-operated benchtop cutaway model for the demonstration of the functionality of a variable speed belt drive
Dimensions and weight:
L x w x h: 350 x 300 x 200 mm
Weight: approx. 7.5 kg.



Laboratory Instrument India