



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cutaway Model Multiple-Disc Clutch

Product Code :
LBNY-0005-10100065



Description :

The mechanical functionality of all components is retained in full. Each of the cutaway models is securely mounted on a base plate to which handles are attached for carrying purposes. The range consists of a series of different cutaway models of drive system components. All components are industrially relevant and have been prepared such that individual active elements are clearly visible to the user. All mechanisms are hand-driven.

Technical Specification :

Industrially-relevant part, fully functional cutaway model

Solid metal base plate, with handles

Hand-operated benchtop cutaway model for the demonstration of the functionality of an multiple-disc clutch

Dimensions and weight:

L x w x h: 350 x 300 x 200 mm

Weight: approx. 6 kg.



Laboratory Instrument India