



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Roller Bearings Kit

**Product Code :**  
LBNY-0005-10100066



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

The kit presents various roller bearings. The bearings are selected for one size of shaft. The kit is clearly laid out on a plastic tray. This instructional kit is used for demonstration and information purposes. The well-structured instructional material enhances the informational value of the kit. No provision is made for conducting exercises or experiments.

**Technical Specification :**

All parts clearly laid out on a plastic tray  
Multiple trays stackable  
Roller bearings demonstration kit  
2 axial bearings / 11 radial bearings  
13 roller bearings, clearly laid out: 5 roller bearings and 8 ball bearings

**Technical data:**

Bearing dimensions  
Outside diameter: d=35, 40, 42, 47, 52mm  
Inside diameter: d=20mm  
Bearing width: h=8, 10, 12, 14, 15, 18, 47mm  
Dimensions and weight:  
L x w x h: 500x350x110mm (tray)  
Weight: approx. 6kg.



**Laboratory Instrument India**