



Email : sales@laboratoryinstrumentindia.com

Product Name :
Roller Bearings Kit

Product Code :
LBNY-0005-10100066



www.laboratoryinstrumentindia.com

Description :

The kit presents various roller bearings. The bearings are selected for one size of shaft. The kit is clearly laid out on a plastic tray. This instructional kit is used for demonstration and information purposes. The well-structured instructional material enhances the informational value of the kit. No provision is made for conducting exercises or experiments.

Technical Specification :

All parts clearly laid out on a plastic tray
Multiple trays stackable
Roller bearings demonstration kit
2 axial bearings / 11 radial bearings
13 roller bearings, clearly laid out: 5 roller bearings and 8 ball bearings

Technical data:

Bearing dimensions
Outside diameter: d=35, 40, 42, 47, 52mm
Inside diameter: d=20mm
Bearing width: h=8, 10, 12, 14, 15, 18, 47mm
Dimensions and weight:
L x w x h: 500x350x110mm (tray)
Weight: approx. 6kg.



Laboratory Instrument India