



Email : sales@laboratoryinstrumentindia.com

Product Name :

Airconditioner (1Hp Split Type, Inverter/ R410a)

Product Code :

LBNY-0002-690001



Description :

Airconditioner (1Hp Split Type, Inverter/ R410a), Airconditioner (1Hp Split Type, Inverter/ R410a) Equipments, Airconditioner (1Hp Split Type, Inverter/ R410a) Tools, Airconditioner (1Hp Split Type, Inverter/ R410a) Tool Kits, Airconditioner (1Hp Split Type, Inverter/ R410a) Manufacturers, Airconditioner (1Hp Split Type, Inverter/ R410a) Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Tunisia, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Mauritania, Eritrea, Gambia, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia

Technical Specification :

The bio-filter has two components in it i.e., Biological enzyme and eco filter.

The eco filter filters the airborne dust particles and the biological enzyme avoids re-pollution and eliminates bacteria by dissolving their cell wall, thus giving you healthy air.

The dehumidification maintains comfortable room temperature and also controls humidity thereby controlling the growth of mold and other micro organisms inside the system.

The self diagnosis feature detects any malfunctioning of the AC and shuts off the AC immediately.

This is equipped with anti-rust cabinet that prevents rust and prolongs the life of the AC and ensures durability.

In Super Turbo mode the AC works faster and cools quicker with minimal pulling of power and is employed for optimum output.

The Sleep mode increases the thermostat temperature without affecting your comfort and will also help save electricity.



Laboratory Instrument India