



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Air Conditioner (2Hp Split Type, Inverter/ R410a)

**Product Code :**  
LBNY-0002-690002



**Description :**

Air Conditioner (2Hp Split Type, Inverter/ R410a)

**Technical Specification :**

Manual Horizontal Airflow Direction Control  
Automatic Operation Mode (Cooling)  
Wide and Long Airflow Vane  
Odour Removing Functions  
12-Hour ON/OFF Setting Timer  
LCD Wireless Remote Controller.

Air Conditioner (2Hp Split Type, Inverter/ R410a), Air Conditioner (2Hp Split Type, Inverter/ R410a) Equipments, Air Conditioner (2Hp Split Type, Inverter/ R410a) Tools, Air Conditioner (2Hp Split Type, Inverter/ R410a) Tool Kits, Air Conditioner (2Hp Split Type, Inverter/ R410a) Manufacturers, Air Conditioner (2Hp Split Type, Inverter/ R410a) Suppliers from India, China, Kenya



**Laboratory Instrument India**

