



Email : sales@laboratoryinstrumentindia.com

Product Name :
Air Conditioner (2Hp Split Type, Inverter/ R410a)

Product Code :
LBNY-0002-690002



Description :

Air Conditioner (2Hp Split Type, Inverter/ R410a)

Technical Specification :

Manual Horizontal Airflow Direction Control
Automatic Operation Mode (Cooling)
Wide and Long Airflow Vane
Odour Removing Functions
12-Hour ON/OFF Setting Timer
LCD Wireless Remote Controller.

Air Conditioner (2Hp Split Type, Inverter/ R410a), Air Conditioner (2Hp Split Type, Inverter/ R410a) Equipments, Air Conditioner (2Hp Split Type, Inverter/ R410a) Tools, Air Conditioner (2Hp Split Type, Inverter/ R410a) Tool Kits, Air Conditioner (2Hp Split Type, Inverter/ R410a) Manufacturers, Air Conditioner (2Hp Split Type, Inverter/ R410a) Suppliers from India, China, Kenya



Laboratory Instrument India

