



Email : sales@laboratoryinstrumentindia.com

Product Name :
Slipping Friction Apparatus

Product Code :
LBNY-0005-1010091



Description :

Frequency of oscillation of different rods supported by the pulleys enables the coefficient of friction to be determined

Set of 5 rods and tubes in mild steel, brass and aluminium alloy.

Self-contained, bench-mounted apparatus, to measure slipping friction in an oscillating system.

Two contra-rotating `V` pulleys of variable centres distance driven by a variable speed motor.

Stop watch provided

Variables to include weight of specimens.

Technical Specification :

Slipping Friction Apparatus



Laboratory Instrument India