



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Unsymmetrical Bending

**Product Code :**  
LBNY-0005-10100096



**Description :**

2 dial gauges with bracket to record the horizontal and vertical deformation of the beam under load  
Storage system to house the components  
Experimental unit for general and unsymmetrical bending of straight beams  
3 beams: i, l and u profiles  
Clamping flange of beam can be clamped in the pillar free to rotate in any direction  
Clamping flange with angle scale to indicate the angular position of the beam  
Eccentricity of load application point adjustable.

**Technical Specification :**

Unsymmetrical Bending



**Laboratory Instrument India**