



Email : sales@laboratoryinstrumentindia.com

Product Name :
Unsymmetrical Bending

Product Code :
LBNY-0005-10100096



Description :

2 dial gauges with bracket to record the horizontal and vertical deformation of the beam under load
Storage system to house the components
Experimental unit for general and unsymmetrical bending of straight beams
3 beams: i, l and u profiles
Clamping flange of beam can be clamped in the pillar free to rotate in any direction
Clamping flange with angle scale to indicate the angular position of the beam
Eccentricity of load application point adjustable.

Technical Specification :

Unsymmetrical Bending



Laboratory Instrument India