



Email : sales@laboratoryinstrumentindia.com

Product Name :
Demonstration Of Euler Buckling

Product Code :
LBNY-0005-10100097



Description :

Loading of test bars with weights
Test bars cannot be overloaded
White backing wall with grid patterning
Demonstration of elastic buckling
4 steel test bars
Representation of all cases of euler buckling
Storage system to house the components

Technical Specification :

Technical Data:

Test bars
Bar cross-section: 0,5x12mm
Material: steel 1.4310 cold-worked
Quantity: 4
Bar length: 180mm
Buckling loads: approx. 2...32N
Weight set
10x 5N
5x 1N.



Laboratory Instrument India