



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Demonstration Of Euler Buckling

**Product Code :**  
LBNY-0005-10100097



**Description :**

Loading of test bars with weights  
Test bars cannot be overloaded  
White backing wall with grid patterning  
Demonstration of elastic buckling  
4 steel test bars  
Representation of all cases of euler buckling  
Storage system to house the components

**Technical Specification :**

**Technical Data:**

Test bars  
Bar cross-section: 0,5x12mm  
Material: steel 1.4310 cold-worked  
Quantity: 4  
Bar length: 180mm  
Buckling loads: approx. 2...32N  
Weight set  
10x 5N  
5x 1N.



**Laboratory Instrument India**