



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Air Conditioner (3Hp Floor Standing Type Non-Inverter/R410a)

**Product Code :**  
LBNY-0002-690004



**Description :**

Air Conditioner (3Hp Floor Standing Type Non-Inverter/R410a)

**Technical Specification :**

Refrigerant: R 410A

Airflow: selectable fan speed, manual vertical swing & automatic horizontal flow

Usage: Industrial Use, Residential Use, Office Use

Compressor: Hermetically shield scroll

Input Power: 3 phase

Type: Floor Standing, Tower.

An environmentally responsible corporate citizen has taken utmost importance.

Even when selecting air conditioners, a new model featuring refrigerant could be the perfect step in promoting your corporate image..

Air Conditioner (3Hp Floor Standing Type Non-Inverter/R410a), Air Conditioner (3Hp Floor Standing Type Non-Inverter/R410a) Equipments, Air Conditioner (3Hp Floor Standing Type Non-Inverter/R410a) Tools, Air Conditioner (3Hp Floor Standing Type Non-Inverter/R410a) Tool Kits, Air Conditioner (3Hp Floor Standing Type Non-Inverter/R410a) Manufacturers, Air Conditioner (3Hp Floor Standing Type Non-Inverter/R410a) Suppliers from India, China, Kenya



**Laboratory Instrument India**