



Email : sales@laboratoryinstrumentindia.com

Product Name :
Deformation Of Frames

Product Code :
LBNY-0005-10100099



Description :

Statically determinate or statically indeterminate bearing support possible
1 long and 1 short clamping pillar
Roller bearing for statically indeterminate support
Investigation of the deformation of steel frames under load
U-shaped and s-shaped frame
Dial gauges record the deformation of the investigated frame under load
Storage system to house the components
Loading of the frame by weights
2 sets of weights with a movable hook to adjust to any load application point
Experimental setup in frame.

Technical Specification :

Technical Data:

Frame made of steel
Cross-section: 20x10mm
Edge length: 600mm
S-shaped: 600x600mm
U-shaped: 600x600mm
Dial gauges
Measuring range: 0...20mm, graduations: 0,01mm
Weights
2x 1N (hanger)

8x 1N
6x 5N.



Laboratory Instrument India