



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Crank Mechanism Kit

**Product Code :**  
LBNY-0005-101000100



**Description :**

Interchangeable bearing bushes permit simulation of bearing clearance  
Springs simulate gas pressure forces  
Can be used together with gear damage accessory set  
Investigation of the vibrations of crank drives  
Crank drive with adjustable stroke  
Stackable storage system to house the components  
Accessory set for machinery diagnostic training system.

**Technical Specification :**

**Technical Data:**  
Stroke: 50 - 75 - 100mm  
Bearing clearance: 0...1mm  
Pressure spring  
Relaxed length: 170mm  
Spring stiffness:  $R=0,55\text{N/mm}$   
Balance mass total  
490g, rated for operation with 50mm stroke.

---

# Laboratory Instrument India

