



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Photoelastic Experiments With Transmission

**Product Code :**  
LBNY-0005-101000101



**Description :**

All filters with 360° angle scale and marking of the main optical axis  
White light generated using a fluorescent tube and two incandescent lamps  
Monochromatic light (colour yellow) generated using a sodium vapour lamp  
Generation of compression or tension forces by means of a threaded spindle  
Representation of mechanical distribution of stress in photoelastic experiments  
2 plane polarisation filters as polariser and analyser  
2 quarter wave filters to generate circular polarised light  
Filters roller bearing mounted and rotating  
Frame cross-arms height-adjustable  
Complete models in polycarbonate (PC) for demonstration purposes available as accessories.

**Technical Specification :**

Photoelastic Experiments With Transmission



**Laboratory Instrument India**

