



Email : sales@laboratoryinstrumentindia.com

Product Name :
Electromechanical Vibrations Kit

Product Code :
LBNY-0005-101000103



Description :

Power display on display and control unit of base
Asynchronous motor with adjustable gap
Asymmetric magnetic field by winding with shut-off facility
Investigation of vibration behaviour of an electric motor
Variable speed via frequency converter of base unit
Speed display on display and control unit of base
Stackable storage system to house the components
Eccentricity of armature adjustable via position of bearing cover
Accessory set for machinery diagnostic training system.

Technical Specification :

Technical Data:
Asynchronous motor with variable speed
Nominal power output: 370W
Eccentricity of armature: 0...0,2mm
Speed range: 100...6000min⁻¹.

Laboratory Instrument India

