



Email : sales@laboratoryinstrumentindia.com

Product Name :
Verification Of Stress Hypotheses

Product Code :
LBNY-0005-101000104



Description :

Technical Data:

Deformation

Measuring range: 0...10mm

Graduations: 0,01mm

Set of weights to place load on specimens

1x 2N (hanger), 1x 1N, 1x 2N, 1x 4N, 2x 8N

Set of weights to compensate for the load and the loading plate

1x1n, 2x2n, 1x4n, 2x8n

Lever arm: 100mm

Specimens

Length: 49mm

Clamping length: 11,5mm

Specimen diameter in measurement

Cross-section: $d=4\text{mm}$.

Technical Specification :

Verification Of Stress Hypotheses

Laboratory Instrument India

