



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Verification Of Stress Hypotheses

**Product Code :**  
LBNY-0005-101000104



**Description :**

**Technical Data:**

Deformation  
Measuring range: 0...10mm  
Graduations: 0,01mm  
Set of weights to place load on specimens  
1x 2N (hanger), 1x 1N, 1x 2N, 1x 4N, 2x 8N  
Set of weights to compensate for the load and the loading plate  
1x1n, 2x2n, 1x4n, 2x8n  
Lever arm: 100mm  
Specimens  
Length: 49mm  
Clamping length: 11,5mm  
Specimen diameter in measurement  
Cross-section:  $d=4\text{mm}$ .

**Technical Specification :**

Verification Of Stress Hypotheses

---

# Laboratory Instrument India

