



Email : sales@laboratoryinstrumentindia.com

Product Name :
Photoelastic Demonstration

Product Code :
LBNY-0005-101000105



Description :

Photoelastic experiments with an overhead polariscope
All filters rotatable in any direction on the horizontal plane
Plane or circular polarised light possible
Green filter for monochromatic light
Polariser and analyser each comprising a polarisation filter and quarter wave filter
Filter enclosed, with stress-free glazing
8 different polycarbonate (pc) models supplied
Load application device with force gauge for compressive and tensile loading
Storage system to house the components.

Technical Specification :

Photoelastic Demonstration



Laboratory Instrument India