



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Photoelastic Demonstration

**Product Code :**  
LBNY-0005-101000105



**Description :**

Photoelastic experiments with an overhead polariscope  
All filters rotatable in any direction on the horizontal plane  
Plane or circular polarised light possible  
Green filter for monochromatic light  
Polariser and analyser each comprising a polarisation filter and quarter wave filter  
Filter enclosed, with stress-free glazing  
8 different polycarbonate (pc) models supplied  
Load application device with force gauge for compressive and tensile loading  
Storage system to house the components.

**Technical Specification :**

Photoelastic Demonstration



**Laboratory Instrument India**