



Email : sales@laboratoryinstrumentindia.com

Product Name :
Buckling Behaviour Of Bars

Product Code :
LBNY-0005-101000107



Description :

Investigation and testing of all relevant buckling cases
Test bars in various lengths and materials
Force measurement using a hydraulic load cell
Measurement of lateral deflection by dial gauge
Test bar ends pinned or fixed
Spindle to apply forces
Verification of the euler theory of buckling
Experiments in horizontal or vertical orientation
Transverse load application device generates shear forces
Storage system to house the components
Additional experiments with supplementary set.

Technical Specification :

Buckling Behaviour Of Bars



Laboratory Instrument India

