



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Buckling Behaviour Of Bars

**Product Code :**  
LBNY-0005-101000107



**Description :**

Investigation and testing of all relevant buckling cases  
Test bars in various lengths and materials  
Force measurement using a hydraulic load cell  
Measurement of lateral deflection by dial gauge  
Test bar ends pinned or fixed  
Spindle to apply forces  
Verification of the euler theory of buckling  
Experiments in horizontal or vertical orientation  
Transverse load application device generates shear forces  
Storage system to house the components  
Additional experiments with supplementary set.

**Technical Specification :**

Buckling Behaviour Of Bars



**Laboratory Instrument India**

