



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Refrigerator (1 door, 92 L, R600a)

**Product Code :**  
LBNY-0002-690005



**Description :**

Refrigerator (1 door, 92 L, R600a)

**Technical Specification :**

Direct Cool Refrigerator: cools efficiently using less electricity.

Capacity 92 Litres: perfect for bachelors, small families and compact spaces like studio apartments, dorms and offices.

Energy Rating: 1 star (Energy consumption: 240 KWH)

Anti-Bacterial gasket prevents build-up of bacteria and fungi inside the refrigerator, maintaining utmost hygiene.

Powerful Compressor for fast icing; thick PUF insulation improves retention of cold air.

Refrigerator (1 door, 92 L, R600a), Refrigerator (1 door, 92 L, R600a) Equipments, Refrigerator (1 door, 92 L, R600a) Tools, Refrigerator (1 door, 92 L, R600a) Tool Kits, Refrigerator (1 door, 92 L, R600a) Manufacturers, Refrigerator (1 door, 92 L, R600a) Suppliers from India, China, Kenya



**Laboratory Instrument India**

