



Email : sales@laboratoryinstrumentindia.com

Product Name :
Device For Shearing Experiment

Product Code :
LBNY-0005-101000113



www.laboratoryinstrumentindia.com

Description :

Test specimens made of electrical grade copper
Shearing anvil and tension lug made of hardened steel
Supplementary experiment for the universal material tester for shear loading bolts
Double shear method avoids misleading bending stresses.

Technical Specification :

Technical Data:
Specimen
Diameter X Length: 6x26mm.



Laboratory Instrument India