



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Device For Shearing Experiment

**Product Code :**  
LBNY-0005-101000113



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Test specimens made of electrical grade copper  
Shearing anvil and tension lug made of hardened steel  
Supplementary experiment for the universal material tester for shear loading bolts  
Double shear method avoids misleading bending stresses.

**Technical Specification :**

**Technical Data:**  
Specimen  
Diameter X Length: 6x26mm.



**Laboratory Instrument India**