



Email : sales@laboratoryinstrumentindia.com

Product Name :
Universal Testing Frame

Product Code :
LBNY-0005-101000115



www.laboratoryinstrumentindia.com

Description :

Frame corners formed by joints rigid to bending, each with 9 high-strength bolts
Installation on four adjustable vibration damping bearings
Large loading frame for experiments on construction components from steelwork and civil engineering using bending, tensile and compression forces
Constructed from U 400 steel profiles
Width between frame levels 635mm
Frame opening w x h 4100x1700mm
Test force in central position max. 300kN
Test force off-centre max. 2x200kN

Technical Specification :

Technical Data:

Frame opening w x h: 4100x1700mm
Free width between pillars: 635mm
Permissible testing forces
Central position: 300kn
Profile used: U 400
Off-centre: 2x 200kn.

Laboratory Instrument India

